

Therapy Guide: Body Contouring



Laser Irradiation System: Megalas-AM-800



Therapy Time (Minutes)	Laser Dose (J/cm ²	m ² in			
	AM-800A600	AM-800A660	AM-800A666		
	20 x 100mW	40 x 100mW	60 x 100mW	Frequency (Time/Week)	
	■ 0.0064 J/cm²-sec	■ 0.0128 J/cm²-sec	■ 0.0192 J/cm²-sec	(Time) Weekly	
	1 irradiator	2 irradiators (to same area)	3 irradiators (to same area)		
40	15.4	30.7	46.1		
30	11.5	23.0	34.6	1 ~ 3	
20	7.7	15.4	23.0		
15	5.8	11.5	17.3		

J/cm²-sec = Unit Energy Density (Height=30 cm, Area=314 cm²)

The laser dose is calculated with the application of the quantity of irradiators to the treatment area.

The therapy time is decided by patient's condition and areas. 30 minutes is recommended.

An exercise program follows the laser therapy is required.

Laser Irradiation System: emLas-900



Thomas	Laser Do in Treatment A Height			
Therapy Time (Minutes)	emLas-900A	emLas-900B	Frequency (Time/Week)	
	Wavelength: 660nm		(, o ,	
	20 x 50mW	20 x 100mW		
	0.0032 J/cm ² -sec	0.064 J/cm ² -sec		
40	7.7	15.4		
30	5.8	11.5	1~3	
20	3.8	7.7		
15	2.9	5.8		

J/cm²-sec = Unit Energy Density (Height=30 cm, Area=314 cm²)

The laser dose is calculated with the application of the quantity of irradiators to the treatment area .

The therapy time is decided by patient's condition and areas. 30 minutes is recommended.

An exercise program follows the laser therapy is required.

Laser Irradiation System Comparison

	Mogalas AM 900	emLas-900	
	Megalas-AM-800	ellicas-900	
ltem			
Quantity of laser Irradiator	1 ~ 3	1	
Laser Quantity per Irradiator	20	20	
Time Control	1 ~ 60 Minutes	1 ~ 60 Minutes	
Frequency Control	0 ~ 1,000 Hz	0 Hz (CW)	
Operation Time Counter (for maintainance)	Irreversible Counter	No	
	405nm/100mW		
	520nm/25mW	NA	
Available Lagar Wayalangth/	635nm/100mW		
Available Laser Wavelength/ Unit Laser Power	660nm/50mW	660nm/50mW	
Offit Laser Fower	660nm/100mW	660nm/100mW	
	780nm/100mW	780nm/100mW	
	808nm/100mW	808nm/100mW	
Net Weight	36 Kg	6 Kg	
Shipping Weight	72 Kg	9 Kg	