







Laser Hood Comparison Table

Manufacturer	Konftec			TTAB	Dermotricos	Luce	Sunetics
Model No.	AM-320A	AM-320B	AM-800A666	Revage670	Tricostym	LDS 100	Model G Laser
Reference Picture							
Laser Device Class	Class IIIa	Class IIIb	Class IIIb	Class IIIa	Class IIIb	Class IIIb	Class IIIa
Laser Diode Temp. Rating	RED: 60°C	RED: 60°C IR: 70°C	RED: 60°C	RED: 40°C	RED: 40°C IR: 40°C	RED: 40°C IR: 40°C	RED: 40°C
Laser Diode Grade	Industrial			Commercial	Commercial	Commercial	Commercial
Laser Diode MTBF (Hours)	> 5,000			NA	NA	NA	5,000
Laser Diode	RED 650 nm						107 x 5 mW
	RED 660 nm	40 x 50 mW	30 x 50 mW	60 x 50 mW		36 x 5 mW	36 x 5 mW
	RED 670 nm				24 x 5 mW		
	IR 780 nm		10 x 90 mW				
Total Laser Power (mW)	RED	2,000	1,500	3,000	120	180	180
	IR		900			945	945
	RED + IR	2,000	2,400	3,000		1,125	1,125
Laser Biostimulation							
Time Counter Equipped	Irreversible Counter			NA	NA	NA	Panel Display
Warranty	2 Years			1 Year	1 Year	1 Year	1 Year
Money Back Guarantee	90 Days			NA	NA	NA	NA
Therapy Time (Min.)							
Androgenetic Alopecia	15	15	15	20	20	20	NA
Alopecia Areata	15	15	15		20	20	
Hair Loss Prevention	15	15	15		20	20	
Before/After Hair Transplantation	15	15	15		15	15	
Session (Time/Week)							
Androgenetic Alopecia	2	2	2	Phase-1: 2/week for 12 weeks.	2	2	NA
Alopecia Areata	2	2	2	Phase-2: 1/ week for 12 weeks.	2	2	
Hair Loss Prevention	1	1	1	Phase-3: 1/ every other week for 12 weeks.	1	1	
Before/After Hair Transplantation	2	2	2	Phase-4: 1/month for 3 months.	2	2	