






Laser Hoods Comparison Table

Manufacturer	Konftec			TTAB	Dermotricos	Luce	Sunetics
Model No.	AM-320A	AM-320B	AM-800A666	Revage670	Tricostym	LDS 100	Model G
Reference Picture							
Laser Class	Class IIIa	Class IIIb	Class IIIb	Class IIIa	Class IIIb	Class IIIb	Class IIIa
Laser Diode Temp. Rating	RED: 60°C	RED: 60°C IR: 70°C	RED: 60°C	RED: 40°C	RED: 40°C IR: 40°C	RED: 40°C IR: 40°C	RED: 40°C
Laser Diode Grade	Industrial Grade			Commercial	Commercial	Commercial	Commercial
Laser Module Grade	Industrial Grade			Commercial	Commercial	Commercial	Commercial
Laser Diode MTBF (Hours)	> 5,000	> 5,000	> 5,000	NA	NA	NA	5,000
Laser Diode Wavelength & Power	650 nm						107 x 5 mW
	660 nm	40 x 50 mW	30 x 50 mW	60 x 50 mW		36 x 5 mW	36 x 5 mW
	670 nm				30 x 4 mW		
	780 nm		10 x 90 mW				
Total Laser Power (mW)	940 nm				189 x 5 mW	189 x 5 mW	
	RED	2,000 mW	1,500 mW	3,000 mW	120 mW	180 mW	535 mW
	INFRARED		900 mW			945 mW	
	RED + IR	2,000 mW	2,400 mW	3,000 mW	120 mW	1,125 mW	535 mW
Time Counter Equipped	☑ Irreversible Counter			☒	☒	☒	☑ Panel Display
Warranty	2 Years			1 Year	1 Year	1 Year	1 Year
Money Back Guarantee	90 Days			NA	NA	NA	NA
Therapy Time (Min.)							
Androgenetic Alopecia	15	15	15	20	20	20	NA
Alopecia Areata	15	15	15		20	20	
Hair Loss Prevention	15	15	15		20	20	
After Hair Transplantation	15	15	15		15	15	
Session (Time / Week)							
Androgenetic Alopecia	2	2	2	Phase-1: 2/week for 12 weeks.	2	2	NA
Alopecia Areata	2	2	2	Phase-2: 1/ week for 12 weeks.	2	2	
Hair Loss Prevention	1	1	1	Phase-3: 1/ every other week for 12 weeks.	1	1	
After Hair Transplantation	2	2	2	Phase-4: 1/month for 3 months.	2	2	